WARNUNG:

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WARNING!

Do not chew your mouthguard.

For a mouthguard to fit, hold and protect, it must sit correctly on your teeth and gums as indicated in the instructions.

If a mouthguard does not fit properly and/or loses its fit by chewing, its protective function and durability will be significantly reduced.

Note: Read the instructions manual thoroughly first, then start with the adjustments step by step.

INSTRUCTIONS FOR USE

You will additionally need: water, a cooking pot, a stopwatch, a spoon, a mirror.

1. Bring a pot half full of water to the boil, remove the boiling water from the heat and leave it to stand for approx. 30 seconds.
Place the mouthguard in the water for **exactly 90 seconds**. Use a stopwatch!
2. After exactly 90 seconds, remove the mouthguard from the boiling water with a spoon. Hold the mouthguard briefly under cold water for one second. Then continue with step 3 IMMEDIATELY.
3. Use the mirror and place the mouthguard in your mouth and insert it into the upper row of teeth. The notch of the mouthguard should be in the middle of the two incisors.
4. When the notch is exactly between the two front teeth, bite lightly on the mouthguard and suck out all the water and air. To do this, press your tongue against the inside of the upper row of teeth. You can also gently help with your fingers and press the outer side of the mouthguard lightly against the teeth or gums.
5. Now leave the mouthguard in the adjusted position for 30 seconds. Then remove the mouthguard from your mouth and hold it under cold running water for 30 seconds.
6. If the mouthguard does not fit well after the first fitting, repeat points 1 to 5.
7. If the mouthguard is too large, you can adjust the ends with a sharp knife or scissors.

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